

Lomas
AMBULANCES

RANGE -
- ROVER

APPROVED BY THE ROVER COMPANY LIMITED



HERBERT LOMAS LIMITED HANDFORTH, WILMSLOW, CHESHIRE, GREAT BRITAIN

Telephone: WILMSLOW 25253 (3 Lines)

Telegrams: *AMBULANCE, HANDFORTH*

Telex: 668913

Lomas' Ambulances known



RECOMMENDED INTERIOR LAYOUT FOR OVERSEAS MARKETS

and used throughout the World



ALTERNATIVE LAYOUTS CAN BE ARRANGED TO CUSTOMER'S
OWN SPECIFICATION

SPECIFICATION

Lomas Coachbuilt Ambulance Body mounted to Range Rover Chassis with 10" Wheelbase Extension

RIGHT OF WAY

Built in illuminated Ambulance Sign

GENERAL

In addition to the 10 inch wheelbase extension, the rear end is modified by bolting and welding a mild steel subframe reinforced and incorporating a hinged rear step which folds flush with body floor.

Front doors are in aluminium, the remainder of the body is a fibre glass construction with appropriate swaging to give extra strength and rigidity. Metal framing is bonded between the inner and outer skins of the fibre glass.

One piece sliding window units in toughened glass in body sides. Two thirds of this glass are opaque or dark glass and the top third—clear glass.

One piece large centre door at rear with fixed window in opaque or dark toughened glass.

PATIENTS' COMPARTMENT

If required Patients Compartment can be separated from driving compartment by means of a bulkhead with sliding glass communication window.

INTERIOR LEFT HAND SIDE (Looking at the Vehicle from the Rear)

Lomas combined/seat and stretcher cot trolley ref. L10 with adjustable head and shoulder movement.

Alternatively a Lomas multi-posture combined seat and stretcher cot trolley can be fitted at extra cost. Ref. L10/5. With either of these trolley stretchers, a back rest with hinged arm rests is fitted and secured permanently to the side of the Ambulance.

INTERIOR RIGHT HAND SIDE (Looking at the Vehicle from the Rear)

The Lomas combined instantaneous seat and stretcher platform is fitted longitudinally and can carry 3 Side Facing sitting case Patients. Alternatively, this Equipment can be quickly converted when required to carry a recumbent Patient by means of a 'wind-out' channel to carry a collapsible stretcher. Stretcher is stowed underneath this seat when not in use.

RECOMMENDED INTERIOR LAYOUT FOR EXPORT

INTERIOR LEFT HAND SIDE (Looking at the Vehicle from the Rear)

LOMAS PATIENT EASY LOADING STRETCHER TRAY fitted complete with automatic locking device. Lightweight alloy stretcher with adjustable head/shoulder movement. Lomas Bed and Pillow.

With this equipment it is not necessary for attendant to enter vehicle when loading or unloading recumbent Patient into Ambulance. Please note when the Ambulance is fitted with left hand drive, this equipment would be fitted on the opposite side.

INTERIOR RIGHT HAND SIDE (Looking at the Vehicle from the Rear)

LOMAS COMBINED INSTANTANEOUS SEAT AND STRETCHER PLATFORM as described above.

Again if the chassis is fitted with left hand drive, this equipment would be fitted on the left hand side of the Ambulance.

SPARE WHEEL

Spare Wheel is located in Driving Compartment behind Passenger/Attendants Seat.

All stretchers are trimmed in a washable material.

The Ambulance is supplied complete with the following accessories: RED or BLUE Emergency Beacon fitted on Roof.

Rubber Bladed Electric Fan.

At the rear on the Right Hand side is an Adjustable Floodlight, to facilitate loading of patients when the Ambulance is being used in darkness.

Syrot Siren is supplied and fitted.

Two Interior Lights fitted in the Patients Compartment are controlled by means of a master switch in the Drivers Cab.

One Thermos Flask and Holder supplied and fitted. Capacity: 568 litre. The complete Ambulance is insulated with fibre glass wool.

Cab Light is fitted in driving compartment.

Fibre Glass First Aid Box is provided complete with contents.

Fire Extinguisher fitted in cab as standard.

Flooring is covered in good quality hard wearing rubber or lino.

Small lockers are fitted in Cab Roof, with access from Ambulance interior.

Blood Transfusion Bottle Holder and Bracket with alternative position. Two exterior mirrors are also supplied and fitted.

When required provision is also made for Oxygen Bottle to be stored in upright position forward of Lomas Patient Easy Loading Stretcher on Right Hand Side.

Interior panels and all trimming in colour to suit Customer's choice.

Exterior of Ambulance painted in a colour to suit Customer's requirements.

Customers transfers affixed to side of Ambulance.

Lettering extra.

Air Conditioning available as optional extra.

INTERIOR DIMENSIONS

Length 7' 2" (1879 mm)

Width 5' 3" (1600 mm)

Height 4' 8" (1422 mm)

SHIPPING DIMENSIONS

Overall Length 16' 2" (4927 mm)

Overall Width 5' 10" (1778 mm)

Overall Height 7' 0" (2133 mm)

Weight 1 ton 19 cwts. 3 qrs. (2019.3 kgs)

The Rover Company Limited is not responsible for the manufacture of approved equipment and/or bodywork, although they have closely examined specification and design.

Matters concerning Sales (including shipping, delivery, etc.), service or Warranty claims are the responsibility of the manufacturer (whose name and address is given on this leaflet) or his agent.

In cases of difficulty concerning approved equipment, the Special Projects Section of the Land-Rover Engineering Department is available to offer advice.

THIS CONVERSION IS SUPPLIED COMPLETE AND DETAILED INFORMATION ON THE SPECIFICATION
SHOULD BE OBTAINED FROM THE MANUFACTURER.